


Montgomery®

Long-Term Cannula

Used to replace either the Short-Term Cannula or a standard tracheostomy tube to establish a long-term tracheostoma.



- Can be trimmed
- Ring washer for safety

27°-angled flange to ensure that the tube rests horizontally when seated in the trachea - nothing protrudes into the trachea.



Dual-purpose groove along the axis of the tube rests in the “6 o’clock” position for orientation and serves as a drainage vent.

24 sizes available



SIZE	Dimensions (mm)		REF for the following stoma lengths (dimension "C") in millimeters							REF	
	A	B	17-25	25-33	33-41	41-49	49-57	57-65	Plug/Ring Set	15 mm Adaptor	
4	9	37	330425	330433	330441	330449	330457	330465	310406	310404	
6	11	37	330625	330633	330641	330649	330657	330665	310606	310604	
8	12	37	330825	330833	330841	330849	330857	330865	310806	310804	
10	13	37	331025	331033	331041	331049	331057	331065	311006	311004	

Sterile, for single patient use. 1 Cannula, 2 Plug/Ring Sets, 1 IFU per carton. Accessories available separately (31xxxx) are supplied non-sterile.

bosmed.com

 neomedic

