



Moore Tracheostomy Tube

To maintain an adequate airway while providing support in the stenotic cervical or thoracic trachea.

- Four components: Silicone tracheostomy tube with obturator and two inner cannulas
- Flexible to conform to individual anatomies while providing patient comfort
- Made of radiopaque silicone to aid in reducing incrustation and clotting
- Standard long length can be trimmed for individual applications
- Extra inner cannulas available as replacement



Moore Tracheostomy Tube

SIZE	REF	Dimensions (mm)		
		A	B	C
6	MTT-006	109	11	6.5
8	MTT-008	109	12	7.5
10	MTT-010	109	13	8.5

Supplied sterile. Includes outer tube, obturator, low-profile inner cannula and inner cannula with 15 mm connector.

Replacement Inner Cannulas

SIZE	REF	Description
6	MTT-006-CN	15 mm Inner Cannula
	MTT-006-LP	Low-Profile Inner Cannula
8	MTT-008-CN	15 mm Inner Cannula
	MTT-008-LP	Low-Profile Inner Cannula
10	MTT-010-CN	15 mm Inner Cannula
	MTT-010-LP	Low-Profile Inner Cannula

Supplied sterile. 2 Cannulas per carton (1 cannula per pouch, 2 pouches per carton).

